



FURSA PADS

*Most Comfortable Pads*

1



Put On (Wear) your Gloves before touching pads

Vaa Gloves zako kabla haujagusu pedi.

2



Put your water on bucket/basin and add five drops of Dettol soap on it.

Weka maji baridi kwenye ndoo au beseni na ongeza matone matano ya Dettol.

3



Soak your pads for three (3) Minutes in cold water mixed Dettol.

Loweke pedi zako kwa dakika tatu(3) kwenye maji baridi yaliyochanganywa na Dettol.

4



Wash Your Pads Well and Remove Every Stain by Soap.

Osha Vizuri Pedi Zako na Ondoa Madoa kwa kutumia sabuni.

5



Use pegs to dry your pads under sunrays to Remove Germs and Bacteria

Tumia vibanio vya nguo kuanika pedi zako kwenye mwanga wa jua ili kuondoa bacteria na vimelea vingine.

6



Clean Hands Before Using Pad.

Nawa Mikono Yako Kabla ya Kutumia Pedi.

7



Stay With Your Pad Atleast Four Hours(4hrs) Depending On Your Menstrual Flow.

Tumia Pedi Yako Kwa Muda Usiozidi Masaa Manne Au Kadri ya Mzunguko Wako.


8



Change Your Pad Every Few Hours or When It is Soaked with Blood.

Badilisha Pedi Yako Mara Tu Inapolowana.

9



Store Your Pads in a Clean and Dry Place for the Next Menstrual Cycle.

Hifadhi Pedi Zako Sehemu Safi na Kavu kwa Matumizi ya Baadae.

10



Reuse Your Pads Almost 24 Times.

Tumia Tena Kwa Mizunguko ya Hedhi Usiopungua 24.

11



Dispose Well After Maximum Use.

Tupa Sehemu Salama Baada ya Matumizi Elekezi.

12



Contact Us for More Information

Kwa Maelezo Zaidi Wasiliana Nasi